



Rotational Viscosity Standards

| Product Code | Type | mPa.s (cP) 20°C | mPa.s (cP) 25°C | Pack Size |
|----------------|----------|-----------------|-----------------|-----------|
| VIS-RT5-600 | RT5 | 5.431 | 4.972 | 600mL |
| VIS-RT10-600 | RT10 | 10.89 | 9.972 | 600mL |
| VIS-RT50-600 | RT50 | 54.43 | 49.53 | 600mL |
| VIS-RT75-600 | RT75 | 84.55 | 77.03 | 600mL |
| VIS-RT100-600 | RT100 | 111.6 | 101.5 | 600mL |
| VIS-RT250-600 | RT250 | 275.7 | 250.7 | 600mL |
| VIS-RT350-600 | RT350 | 386.6 | 351.5 | 600mL |
| VIS-RT500-600 | RT500 | 561.1 | 509.4 | 600mL |
| VIS-RT1K-600 | RT1000 | 1126 | 1018 | 600mL |
| VIS-RT5K-600 | RT5000 | 5492 | 4962 | 600mL |
| VIS-RT10K-600 | RT10000 | 11142 | 10040 | 600mL |
| VIS-RT12K-600 | RT12000 | 14159 | 12769 | 600mL |
| VIS-RT30K-600 | RT30000 | 33168 | 29907 | 600mL |
| VIS-RT60K-600 | RT60000 | 67226 | 60586 | 600mL |
| VIS-RT100K-600 | RT100000 | 110159 | 99565 | 600mL |

All values stated are nominal - 08/12/2008