



## Sucrose Brix / RI Standards

Product Code	Shelf Life	Brix	Certified Value	Index	Pack Size
SS05	6 Weeks	5	5 +/-0.02%	1.34026	15mL
SS075	6 Weeks	7.5	7.5 +/-0.02%	1.34400	15mL
SS10	6 Weeks	10	10 +/-0.02%	1.34782	15mL
SS112	6 Weeks	11.2	11.2 +/-0.02%	1.34968	15mL
SS115	6 Weeks	11.5	11.5 +/-0.02%	1.35015	15mL
SS12	6 Weeks	12	12 +/-0.02%	1.35092	15mL
SS125	6 weeks	12.5	12.5 +/-0.02%	1.35171	15mL
SS15	6 Weeks	15	15 +/-0.02%	1.35568	15mL
SS20	6 Weeks	20	20 +/-0.02%	1.36384	15mL
SS25	6 Weeks	25	25 +/-0.02%	1.37233	15mL
SS30	6 Weeks	30	30 +/-0.02%	1.38115	15mL
SS35	6 Weeks	35	35 +/-0.02%	1.39032	15mL
SS40	6 Weeks	40	40 +/-0.02%	1.39986	15mL
SS45	6 Weeks	45	45 +/-0.02%	1.40978	15mL
SS50	6 Weeks	50	50 +/-0.02%	1.42009	15mL
SS55	6 Weeks	55	55 +/-0.02%	1.43080	15mL
SS60	6 Weeks	60	60 +/-0.02%	1.44193	15mL
SS00	6 Weeks	-	Pure Water	1.33299	15mL

*All values stated are nominal - 08/12/2008*